

MADE IN EUROPE



This product has been developed and manufactured using both raw materials and with production in Europe.

Choose quality, save money, and protect the environment.

(UK) User Manual

Sliding Material Prima

02_20230801_UK



-
-
-
-

Article no.	HMI no.	Name - Size	UDI-DI/GTIN
MC 001-1023	129867	Sliding Material Prima - roll 10 m	5744001973868
MC 001-1027	130155	Sliding Material Prima - roll 20 m	5744001971536
MC 001-1028	129866	Sliding Material Prima - roll 50 m	5744001973899
MC 001-1024	129865	Sliding Material Prima - roll 100 m	5744001972014

Regulation MDR 2017/745 for medical devices class 1



ISO 12 31 03 21

The Original Nordic
Human Handling Concept



Developed & Manufactured by:



Thank you for choosing a product from master care A/S.

We recommend reading the user manual thoroughly before using the product (-s).

The information in this user manual applies to the shown products provided that the products have not been modified or customized and that the products have been used as recommended in this user manual.

If the customer changes, repairs or uses the product in other combinations than recommended by master care A/S the CE-approval no longer applies and the guarantee becomes void.

In case of doubt, contact master care A/S or our distributor in your area.

MDR (EU) 2017/745:

Products in this user manual are CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation MDR (EU) 2017/745 – Class 1.

The intended use of the product/-s:

The product is to be used bilayered to reduce friction for users with decreased function during e.g. transfer, dressing or training. The product can be used to reduce the risk of shear injuries.

Correct use of the product:

Place the Nylon underneath the points/support areas where friction is to be reduced. Always use the product double layered to achieve maximum functionality and reduction of friction. Always remove the product under the user after use to avoid fall injuries and risk of skin overheating.

Warning:

Always do a risk assessment before using the product especially if used e.g., when sitting on the edge of the bed since this poses a risk of the user sliding/falling on to the floor. Always make sure the user has sufficient support underneath their feet when sitting if a nylon

piece is placed underneath their bottom.

Never leave the product on the floor since this poses a risk of falling.

Risk assessment:

Always make a risk assessment and a functional assessment of the user prior to using the product. This way using the product for transferring or moving a user, is ergonomic and legally sound.

Knowledge of the product/transferring:

It is important that the user of the product has been properly instructed in the use of the product and has read the user manual thoroughly prior to using the product.

master care A/S and our distributors provide education in the use of products and on our website all user manuals can be found.

Life span:

The expected life span of products with nylon sliding material with daily wear and tear and washing averagely is 1-3 years. The products life span depends on the use, users load, washing and drying temperatures and how often it is washed and dried.

Make visual inspections frequently and preferably before every use. Products with insufficient friction reduction, defect seams, broken material or other deficiencies must be discarded. This way you will always keep a sound and healthy work environment for both helper and user.

Guarantee:

master care A/S provide a 2-year guarantee covering material defects and manufacturing defects if the product has been used correctly.

Unintended events and injuries when misusing the products:

Unintended events must always be reported to the appropriate National Authority and to master care A/S or the distributor in the current country.

master care A/S and our distributors take no responsibility for personal injuries or product damages if the product(-s) have been used in contravention of the by the manufacturer recommended intended purpose or proper use of the product.

User Manual - Sliding Material Prima

Table of Contents

Turning and repositioning in bed:	page 2
Getting further up in bed, more active user:	page 2
Getting further up in bed, less active user:	page 2
Move the user higher up ontop of the Master Turner:.....	page 3
Transferring from chair and bed:	page 3
Repositioning with Multi Sheet:.....	page 3
Avoid friction and shear against vulnerable skin:	page 4
Technical Specifikations	page 6



User Manual - Sliding Material Prima

Turning and repositioning in bed:



1. Fold the Nylon Master Prima and place the fold under the user, along the spine or under the draw sheet towards the side of the transfer.



2. Bend the user's legs. Place the user's hands on the opposite shoulder or on the pillow. Return the user to supine lying on the Nylon Master Prima and slide/reposition the user to side lying.



3. Remove the Nylon Master Prima by pulling diagonally from the lower part. If you feel a resistance just pull from the opposite lower corner.

Getting further up in bed, more active user:



4. Fold the Nylon Master Prima and place the hand in the fold. From above slide the Nylon under the user's shoulders – first on one side and afterwards on the other side.



5. Ask the user to bend the knees and place feet on the sheet. By pushing the feet down and lifting the pelvis from the sheet the user slides the upper body upwards. If necessary, support the user's feet.



6. Remove the Nylon Master Prima by pulling on the corner of the lower part of the material.

Getting further up in bed, less active user:



7. Turn the user to side lying position. Fold the Nylon Master Prima and place it with the fold under the pelvis. It must cover the pressure points under the head, shoulders and pelvis as the user is turned back to supine lying. If necessary, straighten out the Nylon before the transfer.



8. Ask the user to bend the knees and place the feet on the sheet. Elevate the foot of the bed to use the force of gravity. If necessary, support the user's feet.
Alternatively: Place the Nylon Master Prima under the sheet and use this as a draw sheet during the transfer.



9. Remove the Nylon Master Prima by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

User Manual - Sliding Material Prima

Move the user higher up onto of the Master Turner:



10. Turn the user to side lying position. Fold the Nylon Master Prima and place it with the fold under the pelvis. It must cover the pressure points under the head, shoulders and pelvis as the user is turned back to supine lying.



11. Fix the Master Turner towards the Mattress Cover. Ask user to bend the knees and place the feet on the Mattress Cover for the user to move independently. Alternatively stabilize the user's feet onto the Master Turner/keep the legs straight, elevate the foot end of the bed and push shin/feet to transfer upwards.



12. Remove the Nylon Master Prima by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

Transferring from chair and bed:



13. Ask the user to lean away from the bed. Fold the Nylon Master Prima and slide it under the elevated part of the pelvis.




14. Place the Master Glide sliding board under the folded Nylon Master Prima to build a bridge and make the transferring easier. Remember to make use of the user's resources during the transfer.



15. Remove the sliding board and afterwards remove the Nylon Master Prima by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

Repositioning with Multi Sheet:



16. Turn the user to side lying position using the Multi Sheet if necessary. connected to the ceiling hoist. Fold the Nylon Master and place it with the fold towards the spine. Make sure it covers all the pressure points i.e., pelvis, shoulders, head.



17. Connect the Multi Sheet to the hoist on the opposite side. Elevate and slide the user into position using the Multi Sheet as a draw sheet.




18. Remove the Nylon Master Prima by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

User Manual - Sliding Material Prima

Avoid friction and shear against vulnerable skin:



19. Always remove the Nylon Master Prima from under the user after repositioning.




20. Always make sure you have a firm grip when you work with Nylon Master Prima.



21. Remove the Nylon Master Prima by pulling on the corner of the lower part of the material. This is to avoid unnecessary friction against and tissue injuries to the user's skin.

Product Quality - Functionality - Quality of Life



User Manual - Sliding Material Prima

Technical Specifications

Materials:

Double coated Nylon

ISO certifications:

DS/EN ISO 12952-1:2015

Textiles – Assessing the ignitability of all bedding items - Part 1: Ignition source: smouldering cigarette.

DS/EN ISO 12952-2:2010

Textiles – Assessing the ignitability of all bedding items - Part 2: Ignition source: match-flame equivalent.

Article no. - Product Name:

Size:

Article no. - Product Name:	Size:
MC 001-1023 Sliding Material Prima 10 m	L 10 x W 1,5 m
MC 001-1027 Sliding Material Prima 20 m	L 20 x W 1,5 m
MC 001-1028 Sliding Material Prima 50 m	L 50 x W 1,5 m
MC 001-1024 Sliding Material Prima 100 m	L 100 x W 1,5 m

WARNINGS



Never leave the user on the product when lying on the edge of the bed.



Never leave the product on the floor.



Make sure that all relevant brakes are activated before transferring the user.



CE-certification:

This product is CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation (EU) 2017/745 – class 1.



User Weight:

No maximum weight limit.



Washing:

This product can be washed in a washing machine at max. 85 °C. Gentle washing and dry tumbling significantly prolongs the durability.



Detergents pH level:

This product tolerates detergents / disinfectants with a pH level between 5-9.



Fabric softener:

Never use fabric softener as this will reduce the sliding properties significantly.



Dry Tumbling:

This product may only be dry tumbled gently and at max. 40 °C. Airdrying prolongs the durability significantly.



Surface cleaning:

The surface of this product can be cleaned with a wet cloth and soap.



Disinfectants:

This product tolerates alcoholic disinfectants with a pH level between 5-9.



Ironing:

Do not iron.



Dry cleaning:

Do not dry clean.



Bleach:

Do not bleach.



Maintenance:

Make visual inspections frequently and preferably before every use. Pay attention to torn fabric, seams and straps. Products with high friction or other deficiencies must be discarded.



User Manual:

Always read the User Manual thoroughly before using the product.



Disposal:

Always make sure to follow local guidelines when disposing or discarding the product.

For more information, videos Quick Guide etc., visit www.master-care.dk

